

# Application note: Analysis of BTEX and siloxanes in biogas, syngas




## Context :

Numerous processes have been developed for the valuation of waste and sludge from wastewater treatment plants. The production of methane, in particular from fermentation or pyro-gasification, is a major issue in the constitution of a 100% renewable energy mix in France.

BTEX and siloxanes are often present in biogas before purification because they are widely used in industry and everyday products.

These impurities must be removed prior to any use of the biogas for their deleterious effect on the infrastructure or the health risk they represent.



### Configuration :

Micro GC Fusion, module Rxi-1ms 10m

Injector: Variable volume

Method :

50°C for 40 seconds

Then >> 220 ( 1°C/s) for 40 seconds

Carrier gas: Helium / 25 psi

### Standard :

Trimethylsilanol (TMS) 40 ppm


Benzene 50 ppm

Hexamethyldisiloxane (L2) 45 ppm

Toluene 55 ppm

Xylene 25 ppm

Detection limit = 1 ppm



### Configuration :

Micro GC Fusion, module Stabilwax 10m

Injector : Large volume variable

Method :

80°C for 30 seconds

Then >>100 ( 1°C/s) for 60 seconds

Carrier gas: Helium / 25 psi

### Standard :

Benzene 3 ppm

Toluene 3 ppm

Ethylbenzene 3 ppm

m-xylene 3 ppm

p-xylene 3 ppm

o-xylene 3 ppm